

PACS Migration

Krucom provides PACS migration services for transferring diagnostic images and related data between archives. Our knowledge and experience is a valuable asset through all phases of a successful migration, giving your project a head start.

Assessment


The first stage of a migration project is to assess the systems and requirements involved. Some of the aspects we investigate are:

- What needs to be migrated?
- Which transfer methods are suitable?
- Project timeframe?

Based on the results we present a solution for the migration of e.g. images, overlays, presentation states and reports.

Customisation

Flexibility is an important factor to be able to fulfil individual client needs and requests. Customisation is readily available from Krucom's team of experienced developers. Our standard migrations tools are combined with customised components to efficiently assemble a tailored migration system.



Transfer

The Migration Manager maintains a database of all the data to be moved. It supervises the migration process and registers progress. Customised modules provide the necessary data conversion, RIS synchronisation, prefetch etc.

Verification

When the transfer is complete all anomalies are evaluated and appropriately processed. Finally, the archives are compared to verify that all data is accounted for.

About Krucom

Krucom AB is a system integration company within healthcare. From our office in Lund, Sweden we focus on the European market and cooperate with providers of diagnostic imaging equipment. With over ten years experience and expert knowledge of the DICOM and HL7 standards we offer solutions in the form of products, consultation services and education.